

形成外科用強弯針

1/2 サークル

形成外科
皮膚
一般外科

角針 A
丸針 AR
10本/パック

▶ご注文は
下記の型番で!
(例) A13-70N2 1は白黒糸 2は青糸 3は赤糸
糸の種類 ↑ ↓ 糸の色

針タイプ (原寸)	号数 USP	N ナイロン	B ブレードシルク	T テロンブレード	糸長	型番 (丸針: AR)
A9 AR9	7-0	●	●	●	50cm	A9-70
	6-0	●	●	●		A9-60
	5-0	●				A9-50
	4-0	●				A9-40
A10 AR10	7-0	●	●	●	50cm	A10-70
	6-0	●	●	●		A10-60
	5-0	●				A10-50
	4-0	●				A10-40
A11 AR11	7-0	●	●	●	50cm	A11-70
	6-0	●	●	●		A11-60
	5-0	●	●	●		A11-50
	4-0	●	●	●		A11-40
A12 AR12	7-0	●	●	●	50cm	A12-70
	6-0	●	●	●		A12-60
	5-0	●	●	●		A12-50
	4-0	●	●	●		A12-40
	3-0	●				A12-30

▶ご注文は
下記の型番で!
(例) A13-70N2 1は白黒糸 2は青糸 3は赤糸
糸の種類 ↑ ↓ 糸の色

針タイプ (原寸)	号数 USP	N ナイロン	B ブレードシルク	T テロンブレード	糸長	型番 (丸針: AR)
A13 AR13	7-0	●			50cm	A13-70
	6-0	●	●	●		A13-60
	5-0	●	●	●		A13-50
	4-0	●	●	●		A13-40
	3-0	●				A13-30
A15 AR15	7-0	●			50cm	A15-70
	6-0	●	●	●		A15-60
	5-0	●	●	●		A15-50
	4-0	●	●	●		A15-40
	3-0	●				A15-30
	2-0	●				A15-20
A18 AR18	6-0	●			50cm	A18-60
	5-0	●	●	●		A18-50
	4-0	●	●	●		A18-40
	3-0	●	●	●		A18-30
	2-0	●				A18-20
A19 AR19	5-0	●	●	●	50cm	A19-50
	4-0	●	●	●		A19-40
	3-0	●	●	●		A19-30
	2-0	●				A19-20
A21 AR21	5-0	●			50cm	A21-50
	4-0	●	●	●		A21-40
	3-0	●	●	●		A21-30
	2-0	●				A21-20

形成外科用弱弯針

3/8 サークル

形成外科
皮膚
一般外科

角針 B
丸針 BR

10本/パック

▶ご注文は
下記の型番で！
(例) A13-70N2 1は白黒
2は黒
3は青
糸の種類 ↑ ↓ 糸の色

針タイプ (原寸)	号数 USP	N ナイロン	B ブレードシルク	T テロンブレード	糸長	型番 (丸針: BR)
B9 BR9	7-0	●	●	●	50cm	B9-70
	6-0	●	●	●		B9-60
	5-0	●				B9-50
B10 BR10	7-0	●	●	●	50cm	B10-70
	6-0	●	●	●		B10-60
	5-0	●				B10-50
B11 BR11	7-0	●	●	●	50cm	B11-70
	6-0	●	●	●		B11-60
	5-0	●	●	●		B11-50
	4-0	●				B11-40
B12 BR12	7-0	●	●	●	50cm	B12-70
	6-0	●	●	●		B12-60
	5-0	●	●	●		B12-50
	4-0	●				B12-40
	3-0	●				B12-30
B13 BR13	7-0	●		●	50cm	B13-70
	6-0	●	●	●		B13-60
	5-0	●	●	●		B13-50
	4-0	●	●	●		B13-40
	3-0	●				B13-30
B15 BR15	7-0	●		●	50cm	B15-70
	6-0	●	●	●		B15-60
	5-0	●	●	●		B15-50
	4-0	●	●	●		B15-40
	3-0	●				B15-30
B18 BR18	5-0	●		●	50cm	B18-50
	4-0	●	●	●		B18-40
	3-0	●	●	●		B18-30
	2-0	●				B18-20
B21 BR21	5-0	●		●	50cm	B21-50
	4-0	●	●	●		B21-40
	3-0	●	●	●		B21-30
	2-0	●				B21-20

眼科用弱弯針

3/8 サークル

眼科一般

角針 C
丸針 CR

10本/パック

▶ご注文は
下記の型番で！
(例) A13-70N2 1は白黒
2は黒
3は青
糸の種類 ↑ ↓ 糸の色

針タイプ (原寸)	号数 USP	N ナイロン	B ブレードシルク	SB 軟質ブレード	糸長	型番 (丸針: CR)
C6 CR6	8-0	●	●	●	50cm	C6-80
	7-0	●	●	●		C6-70
	6-0	●		●		C6-60
C7 CR7	8-0	●	●	●	50cm	C7-80
	7-0	●	●	●		C7-70
	6-0	●	●	●		C7-60
	5-0	●		●		C7-50
C9 CR9	8-0	●	●	●	50cm	C9-80
	7-0	●	●	●		C9-70
	6-0	●	●	●		C9-60
	5-0	●		●		C9-50
	4-0	●		●		C9-40
C11 CR11	7-0	●	●	●	50cm	C11-70
	6-0	●	●	●		C11-60
	5-0	●	●	●		C11-50
	4-0	●		●		C11-40
	3-0	●		●		C11-30
C13 CR13	7-0	●	●	●	50cm	C13-70
	6-0	●	●	●		C13-60
	5-0	●	●	●		C13-50
	4-0	●	●	●		C13-40
	3-0	●		●		C13-30
C15 CR15	6-0	●	●	●	50cm	C15-60
	5-0	●	●	●		C15-50
	4-0	●	●	●		C15-40
	3-0	●	●	●		C15-30
	2-0	●		●		C15-20
C18 CR18	5-0	●		●	50cm	C18-50
	4-0	●	●	●		C18-40
	3-0	●	●	●		C18-30
	2-0	●		●		C18-20
C21 CR21	5-0	●		●	50cm	C21-50
	4-0	●	●	●		C21-40
	3-0	●	●	●		C21-30
	2-0	●		●		C21-20

外科用弱弯針

3/8 サークル

外科一般
皮 膚
整 形

角針 E
丸針 ER

10本/パック

▶ご注文は
下記の型番で！
(例) A13-70N2 1は白黒糸 2は黒糸 3は青糸
糸の種類 ↑ ↓ 糸の色

針タイプ (原寸)	号数 USP	N ナイロン	B ブレードシルク	T テロンブレード	糸長	型番 (丸針: ER)
E12 ER12	7-0	●			50cm	E12-70
	6-0	●	●	●		E12-60
	5-0	●	●	●		E12-50
	4-0	●				E12-40
	3-0	●				E12-30
E14 ER14	6-0	●			50cm	E14-60
	5-0	●	●	●		E14-50
	4-0	●	●	●		E14-40
	3-0	●				E14-30
E18 ER18	5-0	●			50cm	E18-50
	4-0	●	●	●		E18-40
	3-0	●	●	●		E18-30
	2-0	●				E18-20
E21 ER21	5-0	●			70cm	E21-50
	4-0	●	●	●		E21-40
	3-0	●	●	●		E21-30
	2-0	●				E21-20
E25 ER25	5-0	●			70cm	E25-50
	4-0	●	●	●		E25-40
	3-0	●	●	●		E25-30
	2-0	●				E25-20

▶ご注文は
下記の型番で！
(例) A13-70N2 1は白黒糸 2は黒糸 3は青糸
糸の種類 ↑ ↓ 糸の色

針タイプ (原寸)	号数 USP	N ナイロン	B ブレードシルク	T テロンブレード	糸長	型番 (丸針: ER)
E29 ER29	4-0	●			70cm	E29-40
	3-0	●	●	●		E29-30
	2-0	●	●	●		E29-20
	1-0	●				E29-10
E33 ER33	4-0	●			70cm	E33-40
	3-0	●	●	●		E33-30
	2-0	●	●	●		E33-20
	1-0	●	●	●		E33-10
E37 ER37	3-0	●	●	●	70cm	E37-30
	2-0	●	●	●		E37-20
	1-0	●	●	●		E37-10
	1	●				E37-1
E41 ER41	2-0	●	●	●	70cm	E41-20
	1-0	●	●	●		E41-10
	1	●				E41-1
	2	●				E41-2
E45 ER45	2-0	●			70cm	E45-20
	1-0	●	●	●		E45-10
	1	●				E45-1
	2	●				E45-2

外科用強弯針

1/2 サークル

外科一般
皮膚
整形

角針 F
丸針 FR
10本/パック

▶ご注文は
下記の型番で！
(例) A13-70N2 1は白黒糸 2は黒糸 3は黒糸
糸の種類 ↑ ↓ 糸の色

針タイプ (原寸)	号数 USP	N ナイロン	B ブレードシルク	T テロンブレード	糸長	型番 (丸針: FR)
F12 FR12	7-0	●			50cm	F12-70
	6-0	●	●	●		F12-60
	5-0	●	●	●		F12-50
	4-0	●				F12-40
	3-0	●				F12-30
F13 FR13	7-0	●			50cm	F13-70
	6-0	●	●	●		F13-60
	5-0	●	●	●		F13-50
	4-0	●	●	●		F13-40
	3-0	●	●	●		F13-30
F14 FR14	7-0	●			50cm	F14-70
	6-0	●	●	●		F14-60
	5-0	●	●	●		F14-50
	4-0	●	●	●		F14-40
	3-0	●	●	●		F14-30
F17 FR17	6-0	●			50cm	F17-60
	5-0	●	●	●		F17-50
	4-0	●	●	●		F17-40
	3-0	●	●	●		F17-30
	2-0	●				F17-20
F21 FR21	5-0	●	●	●	70cm	F21-50
	4-0	●	●	●		F21-40
	3-0	●	●	●		F21-30
	2-0	●	●	●		F21-20
F25 FR25	5-0	●	●	●	70cm	F25-50
	4-0	●	●	●		F25-40
	3-0	●	●	●		F25-30
	2-0	●	●	●		F25-20
	1-0	●				F25-10

▶ご注文は
下記の型番で！
(例) A13-70N2 1は白黒糸 2は黒糸 3は黒糸
糸の種類 ↑ ↓ 糸の色

針タイプ (原寸)	号数 USP	N ナイロン	B ブレードシルク	T テロンブレード	糸長	型番 (丸針: FR)
F29 FR29	4-0	●			70cm	F29-40
	3-0	●	●	●		F29-30
	2-0	●	●	●		F29-20
	1-0	●	●	●		F29-10
F33 FR33	4-0	●			70cm	F33-40
	3-0	●	●	●		F33-30
	2-0	●	●	●		F33-20
	1-0	●	●	●		F33-10
	1	●	●			F33-1
F37 FR37	3-0	●	●	●	70cm	F37-30
	2-0	●	●	●		F37-20
	1-0	●	●	●		F37-10
	1	●	●	●		F37-1
F41 FR41	2-0	●	●	●	70cm	F41-20
	1-0	●	●	●		F41-10
	1	●	●	●		F41-1
	2	●				F41-2
F45 FR45	2-0	●	●	●	70cm	F45-20
	1-0	●	●	●		F45-10
	1	●	●	●		F45-1
	2	●	●	●		F45-2

血管用強弯針

1/2 サークル

血管
心臓血管

丸針 L

10本/パック

▶ご注文は
下記の型番で!
(例) A13-70N2 1は白黒
2は青
3は赤
糸の種類 ↑ ↓ 糸の色

針タイプ (原寸)	号数 USP	N ナイロン	B ブレードシルク	T テロンブレード	糸長	型番
L6	8-0	●	●		50cm	L6-80
	7-0	●	●	●		L6-70
	6-0	●	●	●		L6-60
L8	8-0	●	●		50cm	L8-80
	7-0	●	●	●		L8-70
	6-0	●	●	●		L8-60
	5-0	●	●	●		L8-50
L11	7-0	●	●		50cm	L11-70
	6-0	●	●	●		L11-60
	5-0	●	●	●		L11-50
L14	6-0	●	●	●	50cm	L14-60
	5-0	●	●	●		L14-50
	4-0	●	●	●		L14-40
	3-0	●	●	●		L14-30
L19	5-0	●	●	●	75cm	L19-50
	4-0	●	●	●		L19-40
	3-0	●	●	●		L19-30
L24	5-0	●	●	●	75cm	L24-50
	4-0	●	●	●		L24-40
	3-0	●	●	●		L24-30

血管用弱弯針

3/8 サークル

血管
心臓血管

丸針 M

10本/パック

▶ご注文は
下記の型番で!
(例) A13-70N2 1は白黒
2は青
3は赤
糸の種類 ↑ ↓ 糸の色

針タイプ (原寸)	号数 USP	N ナイロン	B ブレードシルク	T テロンブレード	糸長	型番
M6	8-0	●	●		50cm	M6-80
	7-0	●	●			M6-70
	6-0	●	●	●		M6-60
M8	8-0	●	●		50cm	M8-80
	7-0	●	●			M8-70
	6-0	●	●	●		M8-60
	5-0	●				M8-50
M10	7-0	●	●		50cm	M10-70
	6-0	●	●	●		M10-60
	5-0	●	●	●		M10-50
M12	4-0	●	●		50cm	M10-40
	6-0	●	●	●		M12-60
	5-0	●	●	●		M12-50
	4-0	●	●	●		M12-40
M14	3-0	●			50cm	M12-30
	6-0	●	●	●		M14-60
	5-0	●	●	●		M14-50
	4-0	●	●	●		M14-40
M17	3-0	●	●	●	75cm	M14-30
	5-0	●	●	●		M17-50
	4-0	●	●	●		M17-40
	3-0	●	●	●		M17-30
M20	2-0	●	●	●	75cm	M17-20
	3-0	●	●	●		M20-30
	2-0	●	●	●		M20-20
M23	1-0	●			75cm	M20-10
	3-0	●	●	●		M23-30
	2-0	●	●	●		M23-20
	1-0	●	●	●		M23-10

胃・腸用強弯針







1/2 サークル

胃・腸用
消化器






丸針 N

10本/パック

▶ご注文は
下記の型番で!
(例) A 13-70N2 1は白
2は黒
3は青
糸の種類 ↑ ↓ 糸の色

針タイプ (原寸)	号数 USP	N ナイロン	B ブレッドシルク			糸長	型番
N14 	6-0	●	●			50cm	N14-60
	5-0	●	●				N14-50
	4-0	●	●				N14-40
N16 	5-0	●	●			50cm	N16-50
	4-0	●	●				N16-40
	3-0	●	●				N16-30
N18 	4-0	●	●			75cm	N18-40
	3-0	●	●				N18-30
	2-0	●	●				N18-20
N20 	4-0	●	●			75cm	N20-40
	3-0	●	●				N20-30
	2-0	●	●				N20-20
N23 	4-0	●	●			75cm	N23-40
	3-0	●	●				N23-30
	2-0	●	●				N23-20
N25 	4-0	●	●			75cm	N25-40
	3-0	●	●				N25-30
	2-0	●	●				N25-20

▶ご注文は
下記の型番で!
(例) A 13-70N2 1は白
2は黒
3は青
糸の種類 ↑ ↓ 糸の色

針タイプ (原寸)	号数 USP	N ナイロン	B ブレッドシルク			糸長	型番
N28 	4-0	●	●			75cm	N28-40
	3-0	●	●				N28-30
	2-0	●	●				N28-20
	1-0	●	●				N28-10
N30 	4-0	●	●			75cm	N30-40
	3-0	●	●				N30-30
	2-0	●	●				N30-20
	1-0	●	●				N30-10
N33 	3-0	●	●			75cm	N33-30
	2-0	●	●				N33-20
	1-0	●	●				N33-10
N38 	2-0	●	●			75cm	N38-20
	1-0	●	●				N38-10
	1	●	●				N38-1
N43 	2-0	●	●			75cm	N43-20
	1-0	●	●				N43-10
	1	●	●				N43-1
	2	●	●				N43-2







胃・腸用弱弯針

3/8 サークル

胃・腸用

丸針 NJ

▶ご注文は
下記の型番で!
(例) A13-70N2 1は白
2は黒
3は青
糸の種類 ↑ ↓ 糸の色

針タイプ (原寸)	号数 USP	N ナイロン	B ブレッドシルク	糸長	型番
NJ23 	4-0	●	●	75cm	NJ23-40
	3-0	●	●		NJ23-30
	2-0	●	●		NJ23-20
NJ28 	4-0	●	●	75cm	NJ28-40
	3-0	●	●		NJ28-30
	2-0	●	●		NJ28-20
NJ33 	3-0	●	●	75cm	NJ33-30
	2-0	●	●		NJ33-20
	1-0	●	●		NJ33-10
NJ38 	3-0	●	●	75cm	NJ38-30
	2-0	●	●		NJ38-20
	1-0	●	●		NJ38-10
	1	●	●		NJ38-1

C-45 縫合針消毒容器 円型



縫合針の滅菌・保管用の容器です。
内缶は多孔フタ式のケースで滅菌消毒に使用します。

寸 法: φ60×H20mm
多孔内缶: φ55×H17mm
材 質: 真鍮製 メッキ仕上げ

Handwriting practice area with horizontal dashed lines.