

NATSUME RAT




Together with the Future of Life Science

“The concept of 3Rs^{*1}” proposed by Russell and Burch in 1959 is currently the international principle on performing animal experiments, and this is an essential principle to follow in the life sciences.

^{*1}Replacement, Reduction, and Refinement

Natsume Seisakusho Co.,Ltd. developed and now launch the Natsume Rat, a silicon mannequin for training animal experimental techniques, taking account of how to address the principle of the 3Rs as a manufacturer of animal experiment devices.

We hope to contribute to a reduction in the number of experimental animals used, and Natsume Rat is used by those concerned as much as possible as practical efforts for alternatives to animal experiments in medical, pharmaceutical, and veterinary education and the education for sophisticating animal experiment techniques.



NATSUME SEISAKUSHO CO.,LTD.

Physics and chemistry apparatus / Basic medical apparatus / Pharmaceutical research apparatus / Laboratory animal breeding management apparatus / Medical science apparatus general

Headquarters 2-18-6 Yushima Bunkyo-ku, Tokyo JAPAN
113-8551
Tel: 03-3813-3251 Fax: 03-3815-2002

Osaka Branch 3F Saito Bio Hills Center 7-7-18 Asagi Saito
Ibaraki-shi, Osaka JAPAN 567-0085
Tel: 072-646-9311 Fax: 072-646-9300

<http://www.nazme.co.jp>

NATSUME RATE



KN-590

Rat model for training the animal experimental technique

NATSUME RAT

Specification

Length of body	190mm (7½")
Length of tail	180mm (7")
Sex	Male
Body weight	about 200 g
Material	Silicon


Component parts

Main Body of rat model	× 1
Tail	× 1
Instruction manual	× 1

Related Products (optional)

Simulated blood (50 ml)	× 1
Spare tail	× 1

Structure



Trainable procedures




Retention /
Holding



Oral administration



Endotracheal
intubation



Tail vein injection
Blood collection
through tail vein

- Please, treat Natsume Rat carefully like a living body. Unnatural handling may cause damage. Refer to the instruction manual for the precautions.
- Please, be aware that the appearance and specifications may change for product improvement without a preliminary announcement.