

Individual Identification tools

KN-292 type B , KN-293 type C

Used for mice and rats. A circular truncated-cone shaped probe is designed ear tissue, easily type B and the cone is easily replaceable by screwing.

● KN-292-2-1 type B

Total length: 135mm
Hole diameter: ϕ 1.5mm
Material: SUS

● KN-292-2-2 type B


Total length: 135mm
Hole diameter: ϕ 2.0mm
Material: SUS

● KN-293 type C

For rats and other animals by operating a turntable which has 4 different diameter holes.

Total length: 160mm
Hole diameter: ϕ 1.5, 2.0, 2.5, 3.0mm
Material: SUS

※These ear punches can be autoclaved (high-pressure steam sterilization).



La Pierce

It is the unique ID tools "La Pierce" is available in aluminum or stainless steel. Aluminum one could be chosen from 5 colors (Red, Blue, Light-blue, Green and Black). Stainless steel one is more resistible against scratch than aluminum which is only available in natural material color.

"M series " for Mouse (male tag surface diameter: 3.5mm, needle length: 3.2mm)

M series in Aluminum (A-5056, anodized color, weight 22mg) with female tag

- M-Standard numbering 1 to 20 in 5 colors each (total: 100 tags)
- M-100 numbering 1 to 100 in single color (total: 100 tags)
- M-1000 numbering 1 to 1,000 in single color (total: 1,000 tags)

SM series in Stainless steel (SUS-303, material color, weight 67mg) with female tag.

- SM-100 numbering 1 to 100 (total: 100 tags)
- SM-1000 numbering 1 to 1,000 (total: 1,000 tags)

" R series " for Rat (male tag surface diameter: 3.5mm, needle length: 4.7mm)

R series in Aluminum (A-5056, anodized color, weight 27mg) with female tag

- R-Standard numbering 1 to 20 in 5 colors each (total: 100 tags)
- R-100 numbering 1 to 100 in single color (total: 100 tags)
- R-1000 numbering 1 to 1,000 in single color (total: 1,000 tags)

SR series in Stainless steel (SUS-303, material color, weight: 67mg) with female tag.

- SR-100 numbering 1 to 100 (total: 100 tags)
- SR-1000 numbering 1 to 1,000 (total: 1,000 tags)

Piercing tool (common for "La Pierce" Mouse and Rat)

- PM-3 (Grips with non-slip rubber finish, not applicable for autoclave)
- PM-3S (Grips without non-slip rubber, applicable for autoclave)



PM-3



PM-3S

VI-2208



NATSUME SEISAKUSHO CO.,LTD.

2-18-6, Yushima, Bunkyo-ku, Tokyo JAPAN 113-8551

<http://www.nazme.co.jp> / contact : g-support@nazme.co.jp

