

Individual Identification tools

Ear Punch type B , type C

Truncated cone tip ear punch for mice and rats.
The tip and body are integrated for improved operability.

● KN-292-2-15 type B

Total length: 127mm
Hole diameter: ϕ 1.5mm
Material: SUS

● KN-292-2-20 type B


Total length: 127mm
Hole diameter: ϕ 2.0mm
Material: SUS

● KN-293 type C

For rats and other animals by operating a turntable which has 4 different diameter holes.

Total length: 160mm
Hole diameter: ϕ 1.5, 2.0, 2.5, 3.0mm
Material: SUS

※These ear punches can be autoclaved (high-pressure steam sterilization).



La Pias

Lightweight, easy-to-read tags "La Pias" are ID tools suitable for use with small laboratory animals such as mice and rats. Available in aluminum and stainless steel.
Aluminum : 5 colors (Red, Blue, Light-blue, Green and Black)
Stainless steel : Non-colored (only available in natural material color)
More resistible against scratch compare to aluminum.

"M series " for Mouse (male tag surface diameter: 3.5mm, needle length: 3.2mm)

M series in Aluminum (A-5056, anodized color, weight 22mg) with female tag

- M-Standard numbering 1 to 20 in 5 colors each (total: 100 tags)
- M-100 numbering 1 to 100 in single color (total: 100 tags)
- M-1000 numbering 1 to 1,000 in single color (total: 1,000 tags)

SM series in Stainless steel (SUS-303, material color, weight 67mg) with female tag.

- SM-100 numbering 1 to 100 (total: 100 tags)
- SM-1000 numbering 1 to 1,000 (total: 1,000 tags)

" R series " for Rat (male tag surface diameter: 3.5mm, needle length: 4.7mm)

R series in Aluminum (A-5056, anodized color, weight 27mg) with female tag

- R-Standard numbering 1 to 20 in 5 colors each (total: 100 tags)
- R-100 numbering 1 to 100 in single color (total: 100 tags)
- R-1000 numbering 1 to 1,000 in single color (total: 1,000 tags)

SR series in Stainless steel (SUS-303, material color, weight: 67mg) with female tag.

- SR-100 numbering 1 to 100 (total: 100 tags)
- SR-1000 numbering 1 to 1,000 (total: 1,000 tags)

Piercing tool (common for "La Pias" Mouse and Rat)

- PM-3 (Grips with non-slip rubber finish, not applicable for autoclave)
- PM-3S (Grips without non-slip rubber, applicable for autoclave)



PM-3



PM-3S

V5-2411E



NATSUME SEISAKUSHO CO.,LTD.

2-18-6, Yushima, Bunkyo-ku, Tokyo JAPAN 113-8551

<http://www.nazme.co.jp> / contact : g-support@nazme.co.jp



Our Products

Individual Identification tools

KN-295 Animal Ear Tags

One-touch attachment with dedicated mounting forceps.
Easy-to-read laser marking.

● KN-295-A-1000

Small animal tags No. 1-1000 (Set A)

● KN-295-A-500

Small animal tags No. 1-500 (Set B)

Size : W2.5 X L 6mm

Material: Monel Metal

For Mice, Rats, and Small Animals

● KN-295-B-100

Rabbits and Dogs ear tags No. 1-100

● KN-295-B-1000

Rabbits and Dogs ear tags No. 1-1000

Size : W3.5 X L 15mm

Material: Monel Metal

For Rabbits, dogs and other animals

※ Custom-made numberrings are also available in quantities of 1,000 or more.

Printing: Up to 5 digits (alphabetical characters also available)
(e.g.:A3333)



KN-296 Ear Tag Applicator

Specially designed for small animal ear tags

● KN-296-A

Small animal ear tag applicator for KN-295-A

Total Length: 138mm

Material: SUS



● KN-296-B

Rabbit and dog ear tag applicator for KN-295-B

Total Length: 190mm

Material: SUS with plating finish

