

Individual Identification tools

KN-295 Animal Ear Tags

One-touch attachment with dedicated mounting forceps.
Easy-to-read laser marking.

● KN-295-A-1000

Small animal tags No. 1-1000 (Set A)

● KN-295-A-500

Small animal tags No. 1-500 (Set B)

Size : W2.5 X L 6mm

Material: Monel Metal

For Mice, Rats, and Small Animals

● KN-295-B-100

Rabbits and Dogs ear tags No. 1-100

● KN-295-B-1000

Rabbits and Dogs ear tags No. 1-1000


Size : W3.5 X L 15mm

Material: Monel Metal

For Rabbits, dogs and other animals

※ Custom-made numbrerings are also available in quantities of 1,000 or more.

Printing: Up to 5 digits (alphabetical characters also available)
(e.g.:A3333)



KN-296 Ear tag applicator

Specially designed for small animal ear tags

● KN-296-A

Small animal ear tag applicator for KN-295-A

Total Length: 138mm

Material: SUS



● KN-296-B

Rabbit and dog ear tag applicator for KN-295-B

Total Length: 190mm

Material: SUS with plating finish

