

Animal Marker

Best way to mark the hair



Product page

アニマルマーカ― Animal Marker

Safety-conscious animal markers made by hair dyeing professionals

Uses court dyes specified by the Japanese Ministry of Health, Labor and Welfare.

The product meets the quality standards and testing methods specified by law.

Features

- Markers for easy and safe use on mice and rats.
- Available in 6 colors.
(Blue, Green, Red, Yellow, Orange, Black)
- Strong penetration into hair, and the color does not fade easily by licking.
- About 400 to 600 rats can be processed.
(Depending on the size of the animal and the area to be applied)
- No fading during long-term testing, usually colors remain in good identification for 6~12 weeks.

Animal Maker

Capacity	Units of Sale	Cat No.
10ml	3pcs/box	FG2200
	12pcs/box	F10-
15ml	2pcs/box	FG2000
	12pcs/box	F15-

* When ordering, please include your preferred color code at the end of the cat No..

BLUE : B	GREEN : G	RED : R
YELLOW : Y	ORANGE : D	BLACK : K

Acute oral toxicity studies have been conducted.

The toxicity classification of the animal markers is GHS 5<2000-5000. The LD50 cutoff value was determined to be 2500mg/kg.

(Food and Drug Safety Center)



10ml



15ml

Best way to mark the tail or skin

Oppen
Product page



V4-2410E



Our Products