

Knit gloves with excellent cut resistance

Knit Gloves for small animals

Protect your fingers from Mice, Rats, Hamster, Guinea pigs teeth



Excellent reliability

We use KEVLAR developed by DUPONT in the United States, which has a reputation for excellent cut resistance.

Bite and cut prevention

In order to protect the thumb, index finger, and middle finger that are easily bitten by mice and rats, weaving using KEVLAR prevents injuries and accidents caused by being bitten.



Size / white-S, black-M, purple-L, yellow-LL

Model Number	Size	quantity
FY200DXS	XS	20 pairs/ bag (with aluminum bag)
FY200DS	S	20 pairs/ bag (with aluminum bag)
FY200DM	M	20 pairs/ bag (with aluminum bag)
FY200DL	L	20 pairs/ bag (with aluminum bag)
FY200DLL	LL	20 pairs/ bag (with aluminum bag)
FY500DS	S	10 pairs/ bag (with aluminum bag)
FY500DM	M	10 pairs/ bag (with aluminum bag)
FY500DL	L	10 pairs/ bag (with aluminum bag)
FY500DLL	LL	10 pairs/ bag (with aluminum bag)



The cut resistance of this product means that it is difficult to cut.
It does not mean puncture resistance.

- ① KEVLAR is discolored by sunlight and ultraviolet rays contained in fluorescent lamps.
If the color changes, the strength of the thread may be reduced, so it is recommended to put it in the aluminum bag in which the product is packaged and store it in a dark place after use.
- ② Do not use bleach when washing.
- ③ After washing, please dry in the shade without exposing it to direct sunlight.

V3-2412E



NATSUME SEISAKUSHO CO.,LTD.

2-18-6, Yushima, Bunkyo-ku, Tokyo JAPAN 113-8551
<https://www.nazme.co.jp/> contact : g-support@nazme.co.jp



Our Products

Knit Gloves

Knit Gloves for marmosets

Gloves made of cut-resistant fiber KEVRA to protect fingers that are easily bitten by marmosets. This products is autoclavable.

Standard type FYN504DF-5

Inner gloves FY200D



Long type FYNL504DF-1
FYNL504DF-5



<i>Model Number</i>	<i>Specification</i>	<i>input number</i>
FY200D	Inner gloves S / M / LL	10 pairs/ bag
FYN504DF-5	Standard type Free size	5 pairs/ bag
FYNL504DF-1	Long type 17cm Free size	1 pairs/ bag
FYNL504DF-5	Long type 17cm Free size	5 pairs/ bag

FY500D
Knit Gloves for small animals



Product page

FY
Knit Gloves for marmosets



Product page