


Oral administration needles

KN-349-N Flexible Feeding Needles

A flexible feeding needles that can be easily inserted into mice, rats, guinea pigs, etc. for forced oral administration with extremely low risk of accidents.

Sales unit: 10 pieces/box



Model Number	Tube Dimensions(mm)	SUS Pipe Diameter(mm)
KN-349-M1-N	φ 1.5×L70	outer φ 0.9 x inner φ 0.6
KN-349-RA-N	φ 1.7×L90	outer φ 1.2 x inner φ 0.8
KN-349-RB-N	φ 2.0×L90	outer φ 1.5 x inner φ 1.1
KN-349-RC-N	φ 2.0×L100	outer φ 1.5 x inner φ 1.1
KN-349-RD-N	φ 2.0×L110	outer φ 1.5 x inner φ 1.1



V3-2412E



Product page



NATSUME SEISAKUSHO CO.,LTD.

2-18-6, Yushima, Bunkyo-ku, Tokyo JAPAN 113-8551

<https://www.nazme.co.jp> / contact : g-support@nazme.co.jp



Our Products