

Oral administration needles













KN-348 Oral Dosing Needle

A stainless steel tube with a ball at the end. It can be easily bent by hand.

Material : Needle → Stainless steel

Needle base → Brass chrome plating

Quantity : 6 pieces/bag

Oral administration needles Various												
												
Model number KN-348-	24G-25	24G-38	22G-50	20G-50	20G-70	20G-80	18G-50	18G-80	17G-100	17G-120	15G-80	15G-120
Old model number				KN-348-1	KN-348-2			KN-348-3		KN-384-4		
Gauge No.	24G	24G	22G	20G	20G	20G	18G	18G	17G	17G	15G	15G
shape	Round ball	Round ball	Round ball	Cone	Cone	Cone	Cone	Cone	Cone	Cone	Cone	Cone
Tip size	φ 1.3	φ 1.3	φ 1.3	φ 1.9×4.0	φ 1.9×4.0	φ 1.9×4.0	φ 2.4×4.3	φ 2.4×4.3	φ 3.0×5.0	φ 3.0×5.0	φ 3.0×5.0	φ 3.0×5.0
Needle size (mm)	φ 0.55×25	φ 0.55×38	φ 0.7×50	φ 0.9×50	φ 0.9×70	φ 0.9×80	φ 1.2×50	φ 1.2×80	φ 1.5×100	φ 1.5×120	φ 1.8×80	φ 1.8×120
Approximate weight of (g)	Mouse	~ 14	15 ~ 20	15 ~ 20	20 ~ 25	25 ~ 30	25 ~ 30	25 ~ 30	30 ~ 35	—	—	—
	Rat	—	—	—	—	—	—	50 ~ 70	100 ~ 200	150 ~ 300	150 ~ 300	200 ~ 350
	Hamster	—	—	—	—	—	—	60 ~ 200	—	—	—	—
	Guinea pig	—	—	—	—	—	—	—	—	250 ~ 300	350 ~ 450	250 ~ 300

※ Custom made products to change needle length are also available.

※ The diameter indicates the outer shape.

V1-2208



NATSUME SEISAKUSHO CO.,LTD.

2-18-6, Yushima, Bunkyo-ku, Tokyo JAPAN 113-8551

<http://www.nazme.co.jp> / contact : g-support@nazme.co.jp

