


## Flexible feeding needles

### KN-349 Flexible feeding needles

A flexible feeding needles that can be easily inserted into mice, rats, guinea pigs, etc. for forced oral administration with extremely low risk of accident.

Sales unit: 10 pieces/box



Model Number	Tube Dimensions(mm)	SUS pipe diameter(mm)
KN-349-M1	φ 1.5×L70	Outer φ 0.9 x inner φ 0.6
KN-349-RA	φ 1.7×L90	Outer φ 1.2 x inner φ 0.8
KN-349-RB	φ 2.0×L90	Outer φ 1.5 x inner φ 1.1
KN-349-RC	φ 2.0×L100	Outer φ 1.5 x inner φ 1.1
KN-349-RD	φ 2.0×L110	Outer φ 1.5 x inner φ 1.1

