

# New Product

80th Anniversary



connect to  
the future

May 27  
on sale



600 units  
Limited Edition

NEW

80th Anniversary model

Oppen Orange

KN-297-O-80PJ

Check  
it out!



“For Laboratory Animals Only”  
New color now available for  
individual identification markers!

\* If this product proves popular, we plan to add it to our regular lineup.



<https://www.nazme.co.jp>



80th Anniversary

KN-297-O-80PJ

# Oppen Orange

Laboratory animal's ID marker

Usage Examples

Day1  
ICR



Usage Examples

Day2  
ICR



Usage Examples

Day3  
ICR



Usage Examples

Day7  
ICR



\*The frequency of re-marking varies depending on the strain, individual differences, and age of the animals, etc.  
\*Color appearance may vary depending on the type of lighting.

## [ Features ]

- ◎ The new orange color joins the lineup, bringing the total to six colors.
- ◎ Capable of marking the tail and skin of not only albino mice but also colored mice such as C57BL/6.
- ◎ Animal-safe inks are used. \*Acute oral toxicity testing has been conducted.

Cat No.	Product Name	for Albino Animals	for Colored Animals
KN-297-O-80PJ-6	Oppen 80th Anniversary model Orange [6pcs/pack]	○	○





Laboratory animal's ID marker



# Opppen

KN-297




**Features**

Capable of marking the tail and skin of not only albino mice but also colored mice such as C57BL/6.

Animal-safe inks are used.

**NEW** With its fine tip, even small numbers and characters can be written clearly without smudging, expanding the range of identification options.

**Check !**  
The method of marking the tail with Opppen is suitable for short-term identification. Compared to ear-punching or ear-marking, which involves cutting the animal's ears, this method is less stressful on the animal and easier for the experimenter.



For marking lines on mice, rats, and rabbits. **Opppen 5** colors


Cat No.	Product Name	for Albino Animals	for Colored Animals
KN-297-B-6	Opppen <b>Black</b> [6pcs/pack]	<input type="radio"/>	<input type="radio"/>
KN-297-G-6	Opppen <b>Green</b> [6pcs/pack]	<input type="radio"/>	<input type="radio"/>
KN-297-P-6	Opppen <b>Pink</b> [6pcs/pack]	<input type="radio"/>	<input type="radio"/>
KN-297-LB-6	Opppen <b>Light Blue</b> [6pcs/pack]	<input type="radio"/>	<input type="radio"/>
KN-297-V-6	Opppen <b>Violet</b> [6pcs/pack]	<input type="radio"/>	<input type="radio"/>



For writing numbers and characters on mouse and rat tails **NEW Opppen For Numbers 3** colors

Cat No.	Product Name	for Albino Animals	for Colored Animals
KN-297-N-B-6	Opppen for Numbers <b>Black</b> [6pcs/pack]	<input type="radio"/>	<input type="radio"/>
KN-297-N-P-6	Opppen for Numbers <b>Pink</b> [6pcs/pack]	<input type="radio"/>	<input type="radio"/>
KN-297-N-LB-6	Opppen for Numbers <b>Light Blue</b> [6pcs/pack]	<input type="radio"/>	<input type="radio"/>



Can also be used for marking skin and ears when performing surgery.



\*The frequency of re-marking varies depending on the strain, individual differences, and age of the animals, etc. \*Color appearance may vary depending on the type of lighting.

### Oppen for Numebers (FINE)

3 colors◎Black • Pink • Light Blue



**Final Report**  
Acute oral toxicity test  
of Oppen's ink



**Examples of Use**  
Oppen  
Product page



**Animal Marker**  
Product page

**Suitable for identification  
by staining their hair.**

