



NAZME ID

Laboratory animal's ID marker



Opppen

KN-297



Features

Capable of marking the tail and skin of not only albino mice but also colored mice such as C57BL/6.

Animal-safe inks are used.

NEW With its fine tip, even small numbers and characters can be written clearly without smudging, expanding the range of identification options.

Check !

The method of marking the tail with Opppen is suitable for short-term identification.

Compared to ear-punching or ear-marking, which involves cutting the animal's ears, this method is less stressful on the animal and easier for the experimenter.



Twin marker



THICK

THIN

Single marker



FINE

For marking lines on mice, rats, and rabbits. **Opppen 6 colors**

Cat No.	Product Name	for Albino Animals	for Colored Animals
KN-297-B-6	Opppen Black [6pcs/pack]	○	○
KN-297-G-6	Opppen Green [6pcs/pack]	○	○
KN-297-P-6	Opppen Pink [6pcs/pack]	○	○
KN-297-LB-6	Opppen Light Blue [6pcs/pack]	○	○
KN-297-V-6	Opppen Violet [6pcs/pack]	○	○
KN-297-O-6	Opppen Orange [6pcs/pack]	○	○

For writing numbers and characters on mouse and rat tails

NEW **Opppen** For Numbers **3 colors**

Cat No.	Product Name	for Albino Animals	for Colored Animals
KN-297-N-B-6	Opppen for Numbers Black [6pcs/pack]	○	○
KN-297-N-P-6	Opppen for Numbers Pink [6pcs/pack]	○	○
KN-297-N-LB-6	Opppen for Numbers Light Blue [6pcs/pack]	○	○

Can also be used for marking skin and ears when performing surgery.

ICR

Day 1

Day 2

Day 3

Day 4

Day 7

Black

Green

Pink

Light
Blue

Violet

Orange

* The frequency of re-marking varies depending on the strain, individual differences, and age of the animals, etc. * Color appearance may vary depending on the type of lighting.

Oppen for Numebers (FINE)

3 colors◎Black • Pink • Light Blue



Marking
Sample

THICK
THIN
FINE
Oppen
for Numbers



Final Report

Acute oral toxicity test
of Oppen's ink



Examples of Use

Oppen
Product page



Animal Marker
Product page

**Suitable for identification
by staining their hair.**

【Manufacturer】

V7-2610E



NAZME

Nazme Co. Ltd.
2-18-6, Yushima, Bunkyo-ku, Tokyo 113-8551 JAPAN
Tel : +81-3-3813-3251 Fax : +81-3-3815-2002
HP : <https://www.nazme.co.jp/>



Contact Us