


KN-30651

Laryngoscope I-Shaped Blade for pigs and mini pigs

Specially shaped laryngoscope for pigs.

**A battery inserted in the handle and the tip
can be illuminated with a light bulb.**



株式会社 夏目製作所
<https://www.nazme.co.jp>



Our Products

NATSUME SEISAKUSHO CO.,LTD.

2-18-6, Yushima, Bunkyo-ku, Tokyo JAPAN 113-8551

For Inquiry ☐ g-support@nazme.co.jp

V1-2507E

KN-30651


Laryngoscope I-Shaped Blade for pigs and mini pigs

Specially shaped laryngoscope for pigs.

A battery inserted in the handle and the tip can be illuminated with a light bulb.

【I-Shaped Blade features】

- The I-shaped blade allows visibility from both sides, ensuring a clear view of the inside of the oral cavity.
- The spoon-shaped tip prevents the tongue from falling out.
- The light source has been positioned to allow more light to reach, increasing accuracy and making it brighter and easier to see.
- The laryngoscope I blade can be selected from LED compatible products or light bulb compatible products.



【Specifications】

Material : Brass plated finish

Full length : 260 mm (width of blade tip 20 mm)

Wing length : 145 mm

- KN-30651-L Laryngoscope I-shaped blade for pigs and mini pigs (with LED light bulb)
- KN-30652 Handle for laryngoscope I-shaped blade (KN-30651-LED)
- KN-30651 Laryngoscope I-shaped blade for pigs and mini pigs (with regular light bulb)
- KN-306-D1 Improved Macintosh type blade
Size: L280mm with light bulb
- KN-306-D2 Forager type blade
Size: L300mm with light bulb
- KN-306-D3 Forager type blade
Size: L200mm with light bulb
- KN-306-D4 Handle
- KN-306-D5 Light bulb

